

Synthesis of functionally active proteins containing disulfide bonds using the new PURE system (PUREfrex™)

Takashi Kanamori^{1,2}, Tomoe Murakami², Yuki Hayami¹, Takuya Ueda² (¹GeneFrontier Corp., ²Grad. Sch. of Frontier Sci., Univ. of Tokyo)

PUREfrex™を用いた活性型ジスルフィド結合含有タンパク質の合成

○金森 崇^{1,2}, 村上 朋重², 速水 友紀¹, 上田 卓也² (¹ジーンフロンティア(株), ²東大院・新領域)

Abstract


The PURE system is a reconstituted cell-free protein synthesis system, which consists of only purified factors necessary for transcription, translation and energy regeneration. The PURE system has the unique features that it contains less contaminants such as nucleases and proteases and that composition of the reagents can be easily adjusted in accordance to the purpose. Therefore, the PURE system is now used not only in preparation of the target proteins but also in *in vitro* display technology such as ribosome display (RD). But, we found that the original PURE system yet contained a large quantity of lipopolysaccharide (LPS) from incompletely purified components. We modified the preparation methods of all components that were purified from *E. coli* and developed the new PURE system (PUREfrex™). In PUREfrex, the amount of contaminants such as LPS and RNase is reduced and then both protein synthesis activity and selection efficiency in RD were improved.

However, because PUREfrex includes only translation factors under reduced condition, it is difficult to synthesize proteins containing disulfide bonds in active form. In this study, we examined the composition of reaction mixture and synthesis condition to synthesize functionally active proteins. We used 4 proteins, which form a different number of disulfide bridges, as model proteins. As a result, we succeeded to synthesize all 4 proteins with activities using PUREfrex supplemented with oxidized glutathione (GSSG) and DsbC protein (disulfide bond isomerase) and optimum concentration was different from each target protein.

This result shows that disulfide bonds-containing proteins can be also synthesized in active form using PUREfrex supplemented with optimized concentration of GSSG and DsbC.

1. PURE system

The PURE system is a reconstituted cell-free protein synthesis system, which consists of only purified factors necessary for transcription, translation and energy regeneration.



The original PURE system yet contains contaminants from incompletely purified components.

Advantage

- Low level of contamination
- Easy adjustment of the reagent composition

Application

- High throughput preparation of proteins
- Protein science research
- Translation mechanism
- Folding process
- *in vitro* display
- Ribosome display
- mRNA display

3. Model proteins

- vtPA (truncated tissue plasminogen activator) 9 SS bonds (8 non-consecutive)
 - AppA (*E. coli* acid phosphatase) 5 SS bonds (1 non-consecutive)
 - AP (*E. coli* alkaline phosphatase) 2 SS bonds
 - Fab (Heavy chain + Light chain) 2 SS bonds in each chain
- ### 4. Supplemented factors
- GSSG (oxidized glutathione)
 - DsbC (disulfide bond isomerase in *E. coli*)

2. Development of PUREfrex

Contaminants in PUREfrex are reduced.

	PUREfrex	Original
Proteins		
Tag	None	His-tag
Number of column	3	1
Wash with detergent	+	-
Ribosome		
Wash with detergent	+	-
tRNA		
Wash with detergent	+	-

(Original; Shimizu et al. (2005) Methods, vol.36, p.299)

Table 2-1. Comparison of the preparation methods of components.




Figure 2-1. Comparison of the contaminants in the reagents.

The amount of lipopolysaccharide (LPS), RNase and β -galactosidase (β -Gal) in PUREfrex and the original PURE system was measured. The amount in the original PURE system was set to 100.

Various proteins can be synthesized using PUREfrex.




Figure 2-2. Synthesis of various proteins using PUREfrex.

Various proteins in the figure were synthesized using PUREfrex. The reaction mixture was subjected to SDS-PAGE and the gel was stained with Oriole (Bio-Rad).

PUREfrex is the most suitable system for RD.

	PURE system	<i>E. coli</i>	
	PUREfrex	Original	
Temperature at selection	4-25°C	4-25°C	4°C
RD complex recovery rate	10-30%	2.5%	0.01-0.2%
Enrichment / round	100,000	1,000-10,000	20-250

Table 2-2. Comparison of results of ribosome display among 3 cell-free protein synthesis systems.

5. Synthesis of proteins containing disulfide bonds using PUREfrex SS

Active vtPA can be synthesized in the presence of both GSSG and DsbC.

PUREfrex +GSSG/DsbC
 ↓ + vtPA DNA
 ↓ Incubation at 37°C for 4 h
 Assay of chromogenic substrate cleavage activity
 SDS-PAGE





Figure 5-1. Synthesis of vtPA in the presence of GSSG and DsbC.

vtPA was synthesized using PUREfrex supplemented with the indicated concentration of GSSG and DsbC.
 A, C, 2.5 μ L of the reaction mixture was used in a cleavage assay using chromogenic substrate (SPECTROZYME IPA, American diagnostics). The maximum activity was set to 100.
 B, D, 1 μ L of the reaction mixture was subjected to SDS-PAGE and the gel was stained with SyproOrange (Invitrogen). Synthesized vtPA and supplemented DsbC are indicated by triangles.



Active Fab fragment can be synthesized using PUREfrex with GSSG and DsbC.

PUREfrex +GSSG/DsbC
 ↓ + H chain DNA (+FLAG) / L chain DNA
 ↓ Incubation at 37°C for 4 h
 ↓ to antigen-coated well
 ↓ + anti-FLAG/HRP
 ↓ + TMB




Figure 5-3. Synthesis of Fab fragment using PUREfrex with GSSG and DsbC.

Fab fragment was synthesized using PUREfrex supplemented with 3 mM GSSG and 23 μ g/mL DsbC. The diluted reaction mixture was applied to the antigen-coated well and the bound Fab was detected with anti-FLAG/HRP.

Active phosphatase can be synthesized without DsbC.

PUREfrex +GSSG/DsbC
 ↓ + AppA or AP DNA
 ↓ Incubation at 37°C for 4 h
 Assay of pNPP dephosphorylation activity
 SDS-PAGE




Figure 5-3. Synthesis of acid phosphatase in the presence of GSSG and DsbC.

Acid phosphatase (AppA) was synthesized using PUREfrex supplemented with 3 mM GSSG and the indicated concentration of DsbC.
 A, 2 μ L of the reaction mixture was used in a dephosphorylation assay using p-nitrophenol phosphate (pNPP). The maximum activity was set to 100.
 B, 1 μ L of the reaction mixture was subjected to SDS-PAGE and the gel was stained with SyproOrange. Synthesized AppA and supplemented DsbC are indicated by triangles.




Figure 5-4. Synthesis of alkaline phosphatase.

Alkaline phosphatase (AP) was synthesized using PUREfrex supplemented with the indicated concentration of GSSG in the presence of 94 μ g/mL DsbC or the absence of DsbC.
 A, 2 μ L of the reaction mixture was used in a dephosphorylation assay using pNPP. The maximum activity was set to 100.
 B, 1 μ L of the reaction mixture was subjected to SDS-PAGE and the gel was stained with SyproOrange. Synthesized AP and supplemented DsbC are indicated by triangles.

Figure 5-2. Effect of GSSG and DsbC on the translation of vtPA.

vtPA was synthesized for 2 hours in the presence or absence of GSSG and DsbC. After stopping translation reaction by addition of chloramphenicol (Cm), GSSG and DsbC were added to the indicated samples and all samples were further incubated for 2 hours.
 A, Activity of synthesized vtPA was measured as in Figure 5-1.
 B, Synthesized vtPA was analyzed by SDS-PAGE as in Figure 5-1.

Summary

- We succeeded to synthesize proteins containing disulfide bonds with activities using PUREfrex supplemented with oxidized glutathione and DsbC protein.
- Optimum concentration of oxidized glutathione and DsbC was different from each target protein.

➔ "PUREfrex SS"