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Improvement of translational efficiency by N-terminal codon optimization in the reconstituted cell-free protein synthesis system

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
## 1. PUREfrex<sup>®</sup>; based on the PURE system technology

### Advantage

- Low level of contamination of proteins (including Fab, scFv, protein toxin etc.)
- Easy adjustment of the reagent composition
- PCR products usable as a template DNA

### Application

- High throughput preparation of proteins (including Fab, scFv, protein toxin etc.)
- Protein science research
- In vitro display (Ribosome display, mRNA display etc.)

(Ref: Shimizu *et al.* (2001) *Nat. Biotechnol.*, vol. 19, p. 751.)

## PUREfrex<sup>®</sup> 1.0 PUREfrex<sup>®</sup> 2.0

a regular kit for the synthesis of proteins without disulfide bonds

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a supplement for the synthesis of proteins containing disulfide bonds

### DnaK Mix GroE Mix


a supplement for the synthesis of aggregate-prone proteins

## 2. Experiments


Protein	Reaction mix	Template DNA	Incubation
1. Fab Heavy chain (Herceptin)	PUREfrex <sup>®</sup> 2.0 +DS supplement +DnaK Mix	PCR product	30°C, 4 h
2. PD1 (Programmed Cell Death Protein 1)			
3. PDL1 (Programmed Cell Death Protein 1 Ligand 1)	PUREfrex <sup>®</sup> 2.0 +DS supplement	PCR product	37°C, 4 h
4. GM-CSF (Granulocyte-macrophage colony-stimulating factor)			
5. Citrate Synthase (Yeast Cit1p)			
6. PTGDS (Prostaglandin-H2 D-isomerase)	PUREfrex <sup>®</sup> 2.0	PCR product	37°C, 4 h
7. His-tagged- $\alpha$ -Synuclein			

Protein Synthesis  
↓  
SDS-PAGE (0.5 – 1  $\mu$ L of reaction mixture/lane)  
↓  
Gel stain with Oriole (Bio-Rad) or CBB

## 3. Result 1: N-terminal codon and protein synthesis


## Relationship between protein synthesis and N-terminal sequence


## AT-rich codon vs. GC-rich codon


N-term type	1	2	3	4	5	6	GC(%) 1-6 a.a.
NF001	atg	gaa	ggt	caa	tta	ggt	28%
NF002	atg	gaa	ggt	caa	tta	gta	28%
NF003	atg	gaa	gta	caa	tta	ggt	28%
NF004	atg	gaa	gta	caa	tta	gta	28%
NF377	atg	gag	gtc	cag	ctc	gtc	61%
NF378	atg	gag	gtc	cag	ctc	gtg	61%
NF379	atg	gag	gtc	cag	ctg	gtc	61%
NF380	atg	gag	gtc	cag	ctg	gtg	61%
NF381	atg	gag	gtg	cag	ctc	gtc	61%
NF382	atg	gag	gtg	cag	ctc	gtg	61%
NF383	atg	gag	gtg	cag	ctg	gtc	61%
NF384	atg	gag	gtg	cag	ctg	gtg	61%

## AT-rich codon vs. the most frequently used codon


## 4. Result 2: N-terminal codon optimization of various proteins


## 5. Result 3: Application to Membrane protein synthesis

### N-terminal codon optimization of membrane protein

### Membrane protein synthesis using PUREfrex<sup>®</sup> and Nanodisc

## 6. Conclusion

- N-terminal codon**
- No correlation between codon usage and protein synthesis efficiency.
  - AT-rich codon facilitated the protein synthesis in various proteins including membrane protein.

We recommend the AT-rich codons in N-terminal sequence for high protein expression level.

### On-going

- Investigation of the reason why AT-rich codon of N-terminus facilitates the protein synthesis.

for more information, please contact us.  
URL: www.genefrontier.com  
E-mail: purefrex@genefrontier.com