

# ヒトプロテインジスルフィドイソメラーゼを添加した再構成型無細胞タンパク質合成系(PUREfrex<sup>®</sup>)によるジスルフィド結合タンパク質合成

## Synthesis of proteins containing disulfide bonds using PUREfrex<sup>®</sup> with human protein disulfide isomerase



●松本 令奈、金森 崇 (ジーンフロンティア株式会社)

●Rena Matsumoto and Takashi Kanamori (GeneFrontier Corporation)


### <Abstract>

PURE systemは、大腸菌のタンパク質合成に関与する因子のみから再構成した無細胞タンパク質合成系である。PURE systemにプロテインジスルフィドイソメラーゼ(PDI)として大腸菌DsbCを添加することで、活性を有したジスルフィド結合タンパク質を合成できる。複雑な構造のタンパク質や高等生物のジスルフィド結合タンパク質合成に広く対応するため、本発表ではPUREfrex<sup>®</sup> (改良型PURE system)を用い、DsbCに替えてヒトPDI (hPDI) およびその関連タンパク質を利用できるか検討した。モデル系として、9つのジスルフィド結合を持つ組織プラスミノゲン活性化因子 (tPA) の合成を行った。還元剤に還元型グルタチオン (GSH) を使用したPUREfrexをベースに、酸化型グルタチオン (GSSG) 存在下でhPDIを添加してtPAを合成すると、合成されたtPAはhPDIの添加濃度依存的に活性が増加した。一方、GSSG非存在下でも、hPDIおよびhEro1α (PDI酸化酵素) を添加しtPAを合成した場合、合成されたtPAは同様に高い活性を示した。すなわち、hPDIを機能させるためにはGSSGあるいはhEro1αの添加が必要であることが分かった。また、GSSG存在下でhPDIとDsbCの添加効果を比較したところ、hPDIはDsbCとほぼ同等の効果を示した。以上より、hPDI添加でも、酸化還元環境の調節やhPDI酸化酵素の添加により効率よくジスルフィド結合タンパク質の合成が可能であることが確認された。さらに他の基質タンパク質の合成についても報告する。

### 4. Result 2: Disulfide bonded protein synthesis using hPDI and hEro1α with PUREfrex<sup>®</sup>

#### truncated version of tissue plasminogen activator (tPA)

Organism: *Homo sapiens*  
Synthesized region: 36Ser-40Ile/211Gly-562Pro (+FLAG)  
Length: 368 a.a.  
Molecular weight: 41,072 Da  
No. of disulfide bonds: 9




#### Experimental conditions of protein synthesis


No.	Factors	Basic condition	hPDI	GSSG	hEro1α	DsbC	Template DNA	Temp. & Time
1	hPDI conc.		0-20 μM	2 mM	(-)	(-)	PCR product (1 ng/μL)	30°C, 24 h
2	hPDI + GSSG conc.	PUREfrex <sup>®</sup> 2.1 + 0.5 mM Cys + 4 mM GSH + DnaK mix	10 μM	0-8 mM	(-)	(-)		30°C, 24 h
3	hPDI + hEro1α conc.		10 μM	(-)	0-2 μM	(-)		30°C, 24 h
4	Incubation time		10 μM	2 mM	(-)	(-)		30°C, 0-24 h
5	DsbC		(-)	0 mM/2 mM	(-)	20 μM		30°C, 24 h

●SDS-PAGE (0.25 – 0.5 μL of reaction mix/lane) → Western Blotting (anti-FLAG) or Gel staining (Oriole)  
●tPA activity assay: Assay of chromogenic substrate cleavage activity (2.5 μL PURE reaction mix/ 200 μL assay soln.)


#### 4-1. hPDI concentration





#### 4-2. hPDI + GSSG concentration




#### 4-3. hPDI + hEro1α concentration



#### 4-4. Incubation time



#### 4-5. Comparison of hPDI with DsbC in tPA synthesis



### 5. Result 3: Examples of other disulfide bonded protein synthesis


#### Gaussia Luciferase (GLuc)

Organism: *Gaussia princeps*  
Synthesized region: 181lys-185Asp (+FLAG-His6)  
Length: 187 a.a.  
Molecular weight: 20,407 Da  
No. of disulfide bonds: 5

#### Experimental condition of protein synthesis

No.	Factors	Basic condition	hPDI	GSSG	hEro1α	DsbC	Template DNA	Temp. & Time
1	hPDI		10 μM	(-)	(-)	(-)	PCR product (1 ng/μL)	37°C, 4 h
2	hPDI + hEro1α	PUREfrex <sup>®</sup> 2.1 + 0.5 mM Cys + 4 mM GSH	10 μM	(-)	0.125 μM	(-)		
3	DsbC		20 μM	(-)	(-)	20 μM		
4	GSSG		(-)	1 mM	(-)	(-)		


●SDS-PAGE (0.25 – 0.5 μL of reaction mix/lane) → Western Blotting (anti-FLAG) or Gel staining (Oriole)  
●GLuc activity assay: BioLux Gaussia Luciferase Assay Kit (NEB#E33005) (2 μL PURE reaction mix/50 μL assay soln.)



Combination of hPDI and Ero1α was also effective for active GLuc synthesis. DsbC rather inhibited the activity of synthesized GLuc.

#### Acid phosphatase (AppA)


Organism: *Escherichia coli*  
Synthesized region: 23Gln-432Leu (+FLAG)  
Length: 421 a.a.  
Molecular weight: 45,959 Da  
No. of disulfide bonds: 4



#### Experimental condition of protein synthesis


No.	Factors	Basic condition	hPDI	GSSG	hEro1α	DsbC	Template DNA	Temp. & Time
1	hPDI		2 μM	2 mM	(-)	(-)	PCR product (1 ng/μL)	37°C, 4 h
2	hPDI + GSSG		2 μM	0 / 2 mM	(-)	(-)		
3	hPDI + hEro1α	PUREfrex <sup>®</sup> 2.1 + 0.5 mM Cys + 4 mM GSH	2 μM	(-)	0.2 μM	(-)		
4	DsbC + GSSG		(-)	0 / 2 mM	(-)	4 μM		
5	DnaK mix		0 / 2 μM	0 / 2 mM	0 / 0.2 μM	0 / 4 μM		

●SDS-PAGE (0.5 μL of reaction mix/lane) → Western Blotting (anti-FLAG)  
●AppA activity assay: Assay of PNPP dephosphorylation activity (1 μL PURE reaction mix/100 μL assay soln.)



### 6. Summary

●hPDI can be used as an isomerase with PUREfrex<sup>®</sup>2.1.



PUREfrex<sup>®</sup>2.1 supplemented with hPDI in the presence of GSSG or hPDI oxidase (hEro1α) can synthesize disulfide bonded proteins with their activities. GSSG is unnecessary when oxidizing hPDI is applied with hEro1α.



●Necessary conditions for the synthesis of disulfide bonded proteins with activities are different for each proteins.

Protein name	Basic reaction mix	Reductant	Isomerase	Oxidant/Oxidase
Alkaline Phosphatase (AP)	PUREfrex <sup>®</sup> 2.1 (+ 0.5 mM Cys)	4 mM GSH 2 mM DTT	10 μM hPDI	2 mM GSSG 0.125 μM hEro1α
Truncated version of Tissue Plasminogen Activator (tPA)	PUREfrex <sup>®</sup> 2.1 (+ 0.5 mM Cys) + DnaK mix	4 mM GSH	10 μM hPDI	2 mM GSSG 0.125 μM hEro1α 20 μM DsbC
Gaussia Luciferase (GLuc)	PUREfrex <sup>®</sup> 2.1 (+ 0.5 mM Cys)	4 mM GSH	10 μM hPDI	1 mM GSSG 0.125 μM hEro1α
Acid Phosphatase (AppA)	PUREfrex <sup>®</sup> 2.1 (+ 0.5 mM Cys)	4 mM GSH	2 μM hPDI	2 mM GSSG 0.2 μM hEro1α 4 μM DsbC


In PUREfrex<sup>®</sup>2.1, it is possible to select a suitable isomerase (hPDI etc.) and redox environment depending on the target protein.

For more information, please contact us.  
URL: [www.genefrontier.com](http://www.genefrontier.com)  
E-mail: [purefrex@genefrontier.com](mailto:purefrex@genefrontier.com)

### 1. PUREfrex<sup>®</sup>; based on the PURE system technology



### 2. Oxidative folding of nascent polypeptide in Endoplasmic Reticulum (ER)




#### 2-1. PDIs and related protein used in this study

Protein	UniProt ID	Organism	Length	Molecular weight	PDB ID
hPDI (human protein disulfide isomerase)	P07237	<i>Homo sapiens</i>	492 a.a.	55,424 Da	4EKZ
hEro1α (human endoplasmic oxidoreductin)	Q96HE7	<i>Homo sapiens</i>	446 a.a.	52,122 Da	3AHQ
DsbC (thiol:disulfide interchange protein)	P0AEG6	<i>Escherichia coli</i>	217 a.a.	23,591 Da	1EEJ

### 3. Result 1: Selection of suitable reducing agent

#### Alkaline phosphatase (AP)

Organism: *Escherichia coli*  
Synthesized region: 22Arg-471Lys  
Length: 451 a.a.  
Molecular weight: 47,330 Da  
No. of disulfide bonds: 2



#### Experimental conditions of protein synthesis

No.	Factors	Basic condition	Reductant	hPDI	GSSG	hEro1α	Template DNA	Temp. & Time
1	hPDI		4 mM GSH	10 μM	(-)	(-)	plasmid (10 ng/μL)	37°C, 4 h
2	hPDI + hEro1α		4 mM GSH or 2 mM DTT	10 μM	(-)	0.125 μM		
3	GSSG			(-)	2 mM	(-)		

●SDS-PAGE (0.5 μL of reaction mix/lane) → Gel staining (Oriole)  
●AP activity assay: Assay of PNPP dephosphorylation activity (0.1 μL PURE reaction mix/100 μL assay soln.)

#### 3-1. Difference in synthesized AP activity by reducing agent

